

Colony PCR using Zymolyase

03-07-11, modified by Sigurd Braun (original from Suzanne Komili)

Use:

Diagnostics of yeast (*S. pombe*, *S. cerevisiae*) genomic DNA after chromosomal integration

Requirements:

Zymolyase (e.g. zymolyase 100T, 500 mg, Cat. #120493-1, \$400, **Seikagaki America, Inc.**)

sodium phosphate buffer (0.1 M, pH 7.5), PCR machine

Protocol:

Zymolyase solution:

- Make 2.5 mg/ml zymo solution in 0.1 M sodium phosphate buffer pH 7.5
- Mix by inverting, not all will go into solution
- Spin down and aliquot supernatant (e.g. 100-200 ul) for storage -20° C (alternative: filter sterilize)

note:

in most protocols, the concentration of zymolyase refers to 20T zymolyase; when using 100T zymo, however, less enzyme can be used

Lysis:

- Aliquot 15 ul of zymo solution into PCR tube (can be done in 96well format)
- Scrape small amount of yeast colony (pipette tip size) into PCR tube and swish to resuspend
- Heat to 37° C for 20 min and then 95° C for 5 min in PCR machine
- Mix by inverting, not all will go into solution

PCR:

- Dilute lysis solution 1:10 by adding 150 ul ddH₂O (can now be stored at 4°C or -20°C)
- Use 2.5 ul as template for 25 ul (or 5 in 50 ul) PCR reaction