

Denatured yeast protein extracts by NaOH lysis/TCA precipitation

05-13-06, Sigurd Braun; according to Knop et al. 1999, Yeast 15(10B):963–72

use:

for quick preparation of protein samples (avoid of glass beads and/or phenol and/or boiling)

time: less than 60 minutes

Buffers/Solution:

- **NaOH/BME (Prep under hood fresh each time):** 138.75 μ l 2N NaOH + 11.25 μ l BME per sample
(usually, for 12 samples a mix of 1850 μ l NaOH + 150 μ l BME is prepared!)
- **TCA:** 55% (w/v) trichloroacetic acid in H₂O (stored at 4 deg) – make as follows ;
30ml H₂O + then 55ml OF 100% TCA - + then 15ml H₂O (don't add water to acids – or “don't give acids a drink”)
- **HU buffer:** 200 mM phosphate buffer, pH 6.8, 8 M urea, 5% w/v SDS, 1 mM EDTA, 100 mM DTT, bromophenol blue; stored at –20 deg w/o DTT; **DTT is added freshly from 1 M stock).** (add 80ul of DTT (1M) to 750ul of HU)

cell lysis and protein precipitation:

1. harvest yeast cells corresponding to an OD₆₀₀ = 1 (inoculate from plate with toothpick into 3ml tube - minimal media for S Cere – grow overnight – for S Cere – then measure OD and inoculate into 4ml – (400 x .15/OD) – and grow for 5-6 hours until an OD 0.35-0.50 as measured at a 500ul:500ul dilution= thus 0.6-1 real). You want 1ml with an OD of 1 = so if its below an OD of 1 - then 1000ul/0.6 (for example) OD or above = 1666ul (pipette 833ul twice) – then spin down the pellet and remove supernatant – *You can freeze them at -80 here and resume the next day).*
2. resuspend cell pellet in 1 ml ice-cold water
3. add 150 μ l NaOH/BME to each sample
4. incubate on ice with occasional vortexing for 15 min
5. add 150 μ l TCA to each sample
6. incubate on ice with occasional vortexing for 15 min
7. spin down in a table centrifuge at 4 deg at max. speed (e.g. 20K) for 15 min
8. take off most of supernatant with pipette or with vacuum pump tip/syringe needle
9. spin at max. speed for 5 min, discard rest of SN
10. add 100 μ L ice cooled acetone (flick with finger)
11. spin down at max. speed 15 min., discard SN

12. spin at max. speed for 5 min carefully take off remaining supernatant with precaution (protein pellet becomes slippery after a while)
13. resuspend pellet in 33 μ l HU buffer+DTT (30 μ L HU + 3 μ L DTT), vortex with platform-vortexer 3-5 min
14. optional: when buffer turns yellow due to some residual TCA in the resolved protein pellet, add some (10-20) μ l of 1 M Tris pH 6.8

sample denaturation/SDS PAGE loading:

1. heat denature at 65-70 deg for 10-15 min
(do NOT boil as urea denatures/aggregates at higher temperatures than 70 deg)
2. spin briefly in table centrifuge (RT, 12K, 1 min)
3. load 5-10 μ l of sample (corresponds to 0.1 – 0.2 OD₆₀₀ of yeast cells)