

High-throughput crossing of the *S. pombe* deletion library by Synthetic Genetic Array (SGA) using the Singer RoToR HDA

(updated SGA protocol based on RTR_1v1, by R. R. Barrales as of December 8th 2014)

***S. pombe* library:** The original *Bioneer* deletion library comes in a 96-well format (total 36 plates). For the large-scale crosses, the library has been reassembled into 384 plates (total 9 plates or—after reconfiguration—8 plates, which is the newer version). The library has the mating type *h+*, and the deletions can be selected by Geneticin (G418).

Query strain: Any query strain that is *h-* (or *h90*) and contains another selectable marker (e.g. gene *yΔ::natMX*) can be used for crossing. For selecting haploid spores of only one specific mating type (e.g. *h+*), the Cycloheximide (Cyh) Method by Krogan and co-worker (Roguev et al. 2007) can be applied. In this case, the gene deletion of the query strain needs to be introduced into the original strain (or a derivative) that harbors a *cyh*-resistant allele of the ribosomal RNA gene *rpl42* at its original locus and an additional wild-type *Cyh*-sensitive allele next to the active mating type locus (*h-*). After the cross, haploid spores that contain only the the *Cyh*-resistant allele and are thus *h+* are selected in presence of *Cyh*. For performing silencing assays based on *ura4+* expression/FOA counter-selection, the query strain has additionally a reporter strain with *ura4+* expressed from a heterochromatin locus (i.e. pericentromeric *imr* or *otr*; silent *mat3-EcoRV*). To select for the silent reporter, a selectable marker (usually *hygMX*) is inserted next to the heterochromatin region into a neutral euchromatic locus. In case of the *imrL::ura4+* reporter, *ura4+* is leaky expressed and allows also for selection on EMM plates lacking uracil.

Selection of haploid spores can be done using either using the **Temperature Method** (Dixon et al., PNAS 2008) or the **Cycloheximide Method** (Roguev et al., Nat Methods 2007)

General note for RoToR HDA station settings: On the RoToR touchscreen is a left side tab where several options can be switched on/off (*recycle*, *offset*, *dry mix* or *wet mix*, etc). When loading a profile, these options are not always selected as needed — check first the settings before starting the robot.

Step 1: Generating query strain array and copy of the deletion library

(2 days before mating)

1.1 Create Query Array (384 format): From an over-night culture of the query strain pour 30mL into an empty square plate. Use this plate as source plate for the Singer RoToR HDA station. Several copies (minimum 3 plates) are generated.

*RoToR settings:

- **Select:** Source-Wet 384 bath, Target-Dry 384 plate.
 - Long Pin Pads.
 - Spot Many
- **Load Profile:** Query Array
- **General**
 - **Quantity** Select number of query plates
 - **Pairs** 2 times each (if culture is not saturated, 4 times is recommended)
 - **Recycle** Until repeat (revisit ON)
- **Source**
 - **Pinning**
 - **Speed** 19 mm/sec
 - **Backoff** 0.5 mm
 - **Repeat Pin** 1 time
 - **Wet Mix**
 - **NO WET MIX**
- **Target**
 - **Pinning**
 - **Pinning Pressure** 32%
 - **Speed** 19 mm/sec
 - **Overshoot** 2 mm
 - **Repeat Pin** 1 time
 - **Dry Mix**
 - **NO DRY MIX**

1.2 Create Fresh Library Copy (384 format): replicate all plates of the library (384 format, either assembled in 1-8 or 1-9 plates) onto new YES plates.

*RoToR settings:

▪ Select:	Source-Dry 384 plate, Target-Dry 384 plate.
	<input type="radio"/> Short Pin Pads <input type="radio"/> Replicate Many (in case more than one copy is required)
▪ Load Profile:	SOLID2SOLID SGA
▪ General	
• Quantity	2 (or number of target plates)
• Pairs	1 time each
• Recycle	Until repeat (revisit ON)
▪ Source	
• Pinning	
○ Pinning Pressure	32%
○ Speed	8 mm/sec
○ Overshoot	2 mm
○ Repeat Pin	1 time
• Dry Mix	
○ Clearance	1 mm
○ Diameter	0,5 mm
○ Cycles	1 rotations
▪ Target	
• Pinning	
○ Pinning Pressure	32%
○ Speed	19 mm/sec
○ Overshoot	2 mm
○ Repeat Pin	1 time
• Dry Mix	
○ Clearance	0,5 mm
○ Diameter	1 mm
○ Cycles	1 rotations

Step 2: Mating

2.1 Mix: For mating, replicate each one copy of the query strain and one plate of the deletion library (1, 2...8) onto the *same* SPAS plate (one query strain plate can be used several times, e.g. for up to three library plates).

*RoToR settings:

▪ Select:	Source-Dry 384 plate, Target-Dry 384 plate.
	<input type="radio"/> Short Pin Pads. <input type="radio"/> Mate
▪ Load Profile:	Mating
▪ General	
• Quantity	0
• Pairs	1 time
• Recycle	NO
▪ Source	
• Pinning	
○ Pinning Pressure	32%
○ Speed	19 mm/sec
○ Overshoot	2 mm
○ Repeat Pin	1 time
• Dry Mix	
○ Clearance	0,5 mm
○ Diameter	0,5 mm
○ Cycles	1 rotation

▪	Target		
•	Pinning		
○	Pinning Pressure		32%
○	Speed		19 mm/sec
○	Overshoot		2 mm
○	Repeat Pin	1	time
•	Dry Mix		
○	Clearance		0,5 mm
○	Diameter		1 mm
○	Cycles		1 rotation

2.2 After completing the *Mating* program, adding a **water droplet** on top of cells increases the mating efficiency (Dixon et al. PNAS, 2008). To do this, fill an empty plate with sterile water (20 mL approx.) and run the *H2O Mating* program on the RoToR:

*RoToR settings:

▪	Select:	Source-Wet 384 bath, Target-Dry 384 plate.	
○		Long Pin Pads.	
○		Spot Many	
▪	Load Profile:	H2O Mating	
▪	General		
•	Pairs		1 time
•	Recycle		NO
▪	Source		
•	Pinning		
○	Backoff		0.5 mm
○	Repeat Pin		time
•	Wet Mix		
•	NO WET MIX		
▪	Target		
•	Pinning		
○	Pinning Pressure		32%
○	Speed		19 mm/sec
○	Overshoot		2 mm
○	Repeat Pin		1 time
•	Dry Mix		
○	Clearance		0,5 mm
○	Diameter		0,2 mm
○	Cycles		1 rotation

2.3. Sporulation. Put plates with mated cells in plastic bags and seal them. Incubate **3 to 4 days** at room temperature (25 to 28°C).

Step 3: Haploid Selection

3.1. For **selection of spores** using the **Temperature Method**, put the plates in bags in a **42°C incubator with several containers filled with water (1-2 liters) for 4 days**. Prolonged high-temperature-exposure will kill all cells but spores (i.e. non-mated parental haploid cells, non-sporulated diploid cells).

For selection with the **Cycloheximide (cyh) Method**, continue directly with step 3.2.

3.2. Germination of Spores: For **Temperature Method** replicate spores from SPAS onto YES with 100mg/L of Hyg. For **Cycloheximide (cyh) Method**, replicated spores onto YES with 100 mg/L Cyh). Incubate for 2-3 days at 30°C. During this step, all clones are duplicated changing from 384 to 768 format.

*RoToR settings:

- **Select:** Source-Dry 384 plate Target-Dry 768 plate.
 - Short Pin Pads.
 - 1:2 Array single source.
- **Load Profile:** SOLID2SOLID SGA
- **General**
 - **Quantity** 2 (or number of target plates)
 - **Pairs** 1 time each
 - **Recycle** NO
- **Source**
 - **Pinning**
 - **Pinning Pressure** 32%
 - **Speed** 8 mm/sec
 - **Overshoot** 2 mm
 - **Repeat Pi** 1 time
 - **Dry Mix**
 - **Clearance** 1 mm
 - **Diameter** 0,5 mm
 - **Cycles** 1 rotations
- **Target**
 - **Pinning**
 - **Pinning Pressure** 32%
 - **Speed** 19 mm/sec
 - **Overshoot** 2 mm
 - **Repeat Pin** 1 time
 - **Dry Mix**
 - **Clearance** 0,5 mm
 - **Diameter** 0,3 mm
 - **Cycles** 1 rotations

Step 4: Selection of haploid double mutants

For **Cycloheximide (cyh) Method**, use YES plates supplemented with Cyh (mating type selection), G418 (deletion mutants from library), Hyg (selection marker next to heterochromatin *ura4+* reporter) and another antibiotic (usually Nat for selecting for an additional mutation in the query strain); all antibiotics are at 100 mg/L. For **Temperature Method (cyh) Method**, use YES plates supplemented with G418, Hyg and Nat. **Optional:** When *ura4+* is used as reporter gene within the *imrL* region, an additional EMM-Ura selection step *before* the antibiotics selection step can improve the reproducibility of results.

*RoToR settings:

- **Select:** Source-Dry 768 plate Target-Dry 768 plate.
 - Short Pin Pads.
 - Replicate 768 New2

Note: The option for replicating in the 768 format (*Replicate 768 New2*) is not provided by the RoToR settings but this can be manually programmed and saved to the presets.
- **Load Profile:** SOLID2SOLID 768
- **General**
 - **Quantity** 1
 - **Pairs** 1 time each
 - **Recycle** NO

▪ Source		
• Pinning		
○ Pinning Pressure		10%
○ Speed		19 mm/sec
○ Overshoot		2 mm
○ Repeat Pin		1 time
• Dry Mix		
○ NO DRY MIX		
▪ Target		
• Pinning		
○ Pinning Pressure		10%
○ Speed		19 mm/sec
○ Overshoot		2 mm
○ Repeat Pin		1 time
• Dry Mix		
○ Clearance		0,5 mm
○ Diameter		0,3 mm
○ Cycles		1 rotations

IMPORTANT: If cells have been grown on **EMM-Ura**, it is necessary to use **32% of pinning pressure** and **Dry Mix (0,3 mm diameter)** in *Source* for replica plating to the next selection plate.

Step 5: Silencing assay (with *ura4+* reporter)

Replica-plate from last selection plate onto EMM (control plate), EMM+FOA (1 mg/ml) and EMM-Ura. During this step, the format is changed back from 784 to to 384 format (i.e. splitting each plate with the duplicated colonies from step 3.2. into two plates)

*RoToR settings:

▪ Select:	Source-Dry 768 plate Target-Dry 384 plate.
○	Short Pin Pads.
○	Breakdown
▪ Load Profile:	Silencing Assay From YES
▪ General	
• Quantity	1
• Pairs	1 time each
• Recycle	NO
▪ Source	
• NO Offset	
• Pinning	
○ Pinning Pressure	10%
○ Speed	19 mm/sec
○ Overshoot	2 mm
○ Repeat Pin	1 time
• Dry Mix	
○ NO DRY MIX	
▪ Target	
• Pinning	
○ Pinning Pressure	10%
○ Speed	19 mm/sec
○ Overshoot	2 mm
○ Repeat Pin	1 time
• Dry Mix	
○ NO DRY MIX	

Note: Using these settings, colonies are sometimes not transferred properly from the center of the plate. In this case **increase the pinning pressure to 32 %** and **reduce the spinning speed to 8 mm/sec** for both source and target plate.

Pictures are taking from day 2 to day 4 with a digital camera. Colony sizes are quantified using the software *ht-colony-grid-analyzer-1.1.7* (freeware: <http://ht-col-measurer.sourceforge.net/>)

Required materials:

All quantities (amount of plates) are calculated for one round of large-scale cross (i.e. one query strain with 2 duplicates). Yeast extracts is obtained from Difco, EMM is obtained from ForMedium (EMM without dextrose, #PMD0410). Agar is obtained from Difco (high quality) or Serva (low quality) — we noticed strong variation with respect to growth when using low-quality agar and use for the final steps, in which colony size is quantitatively assessed, only high-quality agar. All antibiotics are used at 100 µg/mL (G418 (geneticin), Invitrogen #1013027; clonNAT (nourseothricin), Werner Bioreagents; hygromycin B, Invitrogen #10687-010; cycloheximide, Singer #C7698). FOA (Thermo Scientific, #R0812) is used at 1 mg/mL.

Media Plates

Plate type	Amount	Step	Comment
YES	3-4	Query strain array (1.1)	Can also be prepared with Nat (or Hyg) to reduce contamination risk
YES	8-10	Library copy (1.2)	Can also be prepared with Gen (G418) to reduce contamination risk
SPAS	8-10	Mating (2.1)	
YES+Hyg	16-20	Germination (3.2)	Temperature method-only
YES+Cyh	16-20	Germination (3.2)	Cycloheximide method-only
EMM-Ura	16-20	Reporter selection (4)	Optional (only for <i>imr::ura4+</i> reporter)
YES+Hyg+Gen+Nat	16-20	Double mutant select (4)	Temperature method-only
YES+Cyh+Hyg+Gen+Nat	16-20	Double mutant select (4)	Cycloheximide method-only
EMM	16-20	Silencing assay (5)	Non-selective conditions
EMM-Ura	16-20	Silencing assay (5)	Ura-positive selection
EMM+FOA	16-20	Silencing assay (5)	Ura-counter selection

Singer RoToR HDA Consumables

dishes and lids (plates): *PlusPlates* #PLU-001

384-pinning pads: *RePads 384 short* #RP-MP-384 (quantity 1000);
RePads 384 long #RP-MP-3L (quantity 200);

Pinning pads	Amount	Step	Comment
384 long	3-4	Query strain array (1.1)	For pinning from liquid cultures, long pins are required
384 short	8	Library copy (1.2)	calculated for 8 plates/library
384 short	16	Mating (2.1)	
384 long	8	Mating (2.2)	Adding water droplet
384 short	16	Germination (3.2)	Temperature method-only
384 short	16	Germination (3.2)	Cycloheximide method-only
384 short	16	Reporter selection (4)	Optional (only for <i>imr::ura4+</i> reporter)
384 short	16	Double mutant select (4)	Temperature method-only
384 short	16	Double mutant select (4)	Cycloheximide method-only
384 short	16	Silencing assay (5)	
384 short	16	Silencing assay (5)	
384 short	16	Silencing assay (5)	
Total number 384 short	120		either temperature or cycloheximide method
Total number 384 long	12		
Total number plates **)	144		

Scheme of screening procedure

